

Diet & Medications

1. It is advised to adhere to usual dietary habits, consuming a variety of foods when possible. Avoid any food that may have potentially resulted in a previous adverse reaction or anaphylaxis. For Celiac Panel, individuals on a gluten-free diet prior to testing may show low or falsely negative serological values.

2. Immunosuppressive drugs such as oral or intranasal corticosteroids (prednisone, beclomethasone, fluticasone, triamcinolone), and topical cortisone suspensions and creams, **may** affect test results for antibody assessment.



16020 Linden Ave N
Shoreline, WA 98133-5672, USA
Phone: 877.318.8728
Phone: 206.365.1256
Fax: 206.363.8790
www.USBioTek.com
Email: cservice@usbiotek.com

To order supplies, call US BioTek Laboratories during our normal business hours: 8:30 a.m. to 5:30 p.m. PT.

Diet & Medications

1. It is advised to adhere to usual dietary habits, consuming a variety of foods when possible. Avoid any food that may have potentially resulted in a previous adverse reaction or anaphylaxis. For Celiac Panel, individuals on a gluten-free diet prior to testing may show low or falsely negative serological values.

2. Immunosuppressive drugs such as oral or intranasal corticosteroids (prednisone, beclomethasone, fluticasone, triamcinolone), and topical cortisone suspensions and creams, **may** affect test results for antibody assessment.



16020 Linden Ave N
Shoreline, WA 98133-5672, USA
Phone: 877.318.8728
Phone: 206.365.1256
Fax: 206.363.8790
www.USBioTek.com
Email: cservice@usbiotek.com

To order supplies, call US BioTek Laboratories during our normal business hours: 8:30 a.m. to 5:30 p.m. PT.

US BioTek Laboratories

Serum Specimen Collection & Shipping Instructions

*Food, Inhalant & Mold Antibody Assessment Panels
Candida Antibodies & Antigen Panel
Celiac Antibody Panel*

Collection Kit Contents

- 2 10-ml Red Top Vacutainer
- 1 10-ml Serum Transfer Tube
- 1 Disposable Transfer Pipette
- 1 Resealable Bag
- 1 Absorbent Pad
- 1 Polyfoam Packer
- 1 Protective Cardboard Box
- 1 Test Requisition Form
- 1 FedEx Express Clinical Pak & Billable Stamp



*Celiac & Specific-IgA Panels are not suitable for breast-fed infants.
Total IgE is not available for age less than 10 years.*

March 2018

US BioTek Laboratories

Serum Specimen Collection & Shipping Instructions

*Food, Inhalant & Mold Antibody Assessment Panels
Candida Antibodies & Antigen Panel
Celiac Antibody Panel*

Collection Kit Contents

- 2 10-ml Red Top Vacutainer
- 1 10-ml Serum Transfer Tube
- 1 Disposable Transfer Pipette
- 1 Resealable Bag
- 1 Absorbent Pad
- 1 Polyfoam Packer
- 1 Protective Cardboard Box
- 1 Test Requisition Form
- 1 FedEx Express Clinical Pak & Billable Stamp



*Celiac & Specific-IgA Panels are not suitable for breast-fed infants.
Total IgE is not available for age less than 10 years.*

March 2018

Important Shipping Instructions

- Clearly print Patient Name (First & Last), Date of Birth & Specimen Collection Date on Serum Transfer Tube and Test Requisition Form.
- Test Requisition Form must also include:
 - Practitioner's Signature · Specimen Storage
 - Requested Test · Payment Information.
- **Do not initiate or request Saturday pickup or delivery from FedEx. Additional shipping charges will be incurred.**
- Missing or incomplete test request information will delay processing and extend turn-around time.
- Complete FedEx Express Billable Stamp and stick onto Clinical Pak.
- Stamps are prepaid for FedEx 2Day delivery for lower 48 states. Contact FedEx for pickup: 1.800.463.3339.
- For clients outside USA & Canada : prepaid shipping does not apply. Client is responsible for the shipping cost.
- For tracking purposes, keep a copy of your tracking number.
- Specimen may be safely shipped at ambient temperature.

Food, Inhalant, Mold & Candida Panels - the lab must receive specimen within 14 days at ambient temperature from date of collection.

Celiac Antibody Panel - the lab must receive specimen within 10 days at ambient temperature from date of collection.

US BioTek Laboratories is NOT a Medicare or Medicaid provider. Our test procedures are not reimbursable by either.

Avoid specimen rejection: label Serum Transfer Tube with Patient Full Name, Date of Birth AND Specimen Collection Date. A specimen tube missing this information will be rejected.

Patient Name on Serum Transfer Tube must exactly match Patient Name on Test Requisition Form. Do not use nicknames or abbreviations.

For IgE Inhalant Panel, patient's zip code is required.

Important Shipping Instructions

- Clearly print Patient Name (first & last), Date of Birth & Specimen Collection Date on Serum Transfer Tube and Test Requisition Form.
- Test Requisition Form must also include:
 - Practitioner's Signature · Specimen Storage
 - Requested Test · Payment Information.
- **Do not initiate or request Saturday pickup or delivery from FedEx. Additional shipping charges will be incurred.**
- Missing or incomplete test request information will delay processing and extend turn-around time.
- Complete FedEx Express Billable Stamp and stick onto Clinical Pak.
- Stamps are prepaid for FedEx 2Day delivery for lower 48 states. Contact FedEx for pickup: 1.800.463.3339.
- For clients outside USA & Canada : prepaid shipping does not apply. Client is responsible for the shipping cost.
- For tracking purposes, keep a copy of your tracking number.
- Specimen may be safely shipped at ambient temperature.

Food, Inhalant, Mold & Candida Panels - the lab must receive specimen within 14 days at ambient temperature from date of collection.

Celiac Antibody Panel - the lab must receive specimen within 10 days at ambient temperature from date of collection.

US BioTek Laboratories is NOT a Medicare or Medicaid provider. Our test procedures are not reimbursable by either.

Avoid specimen rejection: label Serum Transfer Tube with Patient Full Name, Date of Birth AND Specimen Collection Date. A specimen tube missing this information will be rejected.

Patient Name on Serum Transfer Tube must exactly match Patient Name on Test Requisition Form. Do not use nicknames or abbreviations.

For IgE Inhalant Panel, patient's zip code is required.

Specimen Collection Instructions

Serum Requirements:

- Serum volumes are listed on the Test Requisition Form. An insufficient sample may result in the inability to complete testing for all panels ordered. Notification of the panels that could not be completed is issued with the partial items. Follow collection instructions carefully.

There is a \$15 fee for submittal of an unacceptable specimen.

1. Draw the blood utilizing the **Red Top Vacutainer** tubes provided. 1 fully drawn **Red Top Vacutainer** OR SST tube will usually yield about 3-4 ml of serum.
2. Allow the blood to clot at ambient temperature for 30 minutes.
3. Centrifuge at 3000 RPM for 15 minutes.
4. With the **Disposable Transfer Pipette**, transfer serum from the **Red Top Vacutainer** into **Serum Transfer Tube**. Do not transfer red blood cells to the **Serum Transfer Tube**.
5. Secure tightly.
6. Wrap **Serum Transfer Tube** with **Absorbent Pad** and place it into the **Resealable Bag** provided. Seal contents into **Polyfoam Packer** and then into the **Protective Cardboard Box**.
7. Include **Test Requisition Form** inside the **Protective Cardboard Box**.
8. Put the **Protective Cardboard Box** into the **FedEx Express Clinical Pak**. Seal securely.

Hemolyzed, lipemic or icteric sera are rejected.

Specimen Collection Instructions

Serum Requirements:

- Serum volumes are listed on the Test Requisition Form. An insufficient sample may result in the inability to complete testing for all panels ordered. Notification of the panels that could not be completed is issued with the partial items. Follow collection instructions carefully.

There is a \$15 fee for submittal of an unacceptable specimen.

1. Draw the blood utilizing the **Red Top Vacutainer** tubes provided. 1 fully drawn Red Top Vacutainer OR SST tube will usually yield about 3-4 ml of serum.
2. Allow the blood to clot at ambient temperature for 30 minutes.
3. Centrifuge at 3000 RPM for 15 minutes.
4. With the **Disposable Transfer Pipette**, transfer serum from the **Red Top Vacutainer** into **Serum Transfer Tube**. Do not transfer red blood cells to the **Serum Transfer Tube**.
5. Secure tightly.
6. Wrap **Serum Transfer Tube** with **Absorbent Pad** and place it into the **Replaceable Bag** provided. Seal contents into **Polyfoam Packer** and then into the **Protective Cardboard Box**.
7. Include **Test Requisition Form** inside the **Protective Cardboard Box**.
8. Put the **Protective Cardboard Box** into the **FedEx Express Clinical Pak**. Seal securely.

Hemolyzed, lipemic or icteric sera are rejected.